

Anti-PDCD10 antibody (STJ11107392)

STJ11107392

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	11235
Gene Symbol	PDCD10
Uniprot ID	PDC10_HUMAN
Immunogen	
Immunogen Region	
Specificity	A synthetic peptide of human PDCD10
Immunogen Sequence	