

Anti-KCNA1 antibody (500-600 aa) (STJ11107393)

STJ11107393

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	3736
Gene Symbol	KCNA1
Uniprot ID	KCNA1_HUMAN
Immunogen	500-600 aa
Immunogen Region	500-600 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 500-600 of human KCNA1 (Q09470).
Immunogen Sequence	