

Anti-AKAP8L antibody (1-240 aa) (STJ11107394)

STJ11107394

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	26993
Gene Symbol	AKAP8L
Uniprot ID	AKP8L_HUMAN
Immunogen	
Immunogen Region	1-240 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human AKAP8L (NP_055186.3).
Immunogen Sequence	MSYTGFFVQGSETTLQSTYSD TSAQPTCDYGYGTWNSGTNR GYEGYGYGYGYGQDNTTNYG YGMATSHSWEMPSSDTNANT SASGSASADSVLSRINQRLD MVPHLETMMQGGVYGSGGE RYDSYESCDRAVLSELDLY RSGYDYSELDPEMEMAYEGQ YDAYRDQFRMRGNDTFGPRA QGWARDARSGRPMASGYGRM WEDPMGARGQCMSGASRLPS LFSQNIPEYGMFQGMRGG