

Anti-TET3 antibody (STJ11107395)
STJ11107395**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	
Applications	WB/IP
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:200-1:500 IP:1:50-1:100
Formulation	PBS with 0.02% Sodium Azide, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	200424
Gene Symbol	TET3
Uniprot ID	TET3_HUMAN
Immunogen	
Immunogen Region	
Specificity	A synthetic peptide of human TET3
Immunogen Sequence	