

Anti-RPL15 antibody (1-204 aa) (STJ11107398)

STJ11107398

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:200-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	6138
Gene Symbol	RPL15
Uniprot ID	RL15_HUMAN
Immunogen	
Immunogen Region	1-204 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-204 of human RPL15 (NP_002939.2).
Immunogen Sequence	MGAYKYIQELWRKKQSDVMR FLLRVRCWQYRQLSALHRAP RPTRPDKARRLGYKAKQGYV IYRIRVRGGGRKRPVPGAT YGKPVHHGVNQLKFARSLQS VAEERAGRHCALRVLNSYW VGEDSTYKFFEVLIDPFHK AIRRNPDTQWITKPVHKHRE MRGLTSAGRKSRGLGKGHKF HHTIGGSRAAWRRRNTLQL HRYR