

Anti-FGL1 antibody (1-312 aa) (STJ11107401)

STJ11107401

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	2267
Gene Symbol	FGL1
Uniprot ID	FGL1_HUMAN
Immunogen	
Immunogen Region	1-312 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-312 of human FGL1 (NP_004458.3).
Immunogen Sequence	MAKVFSFILVTTALTMGREI SALEDCAQEQMRLRAQVRLLETRVKQQQVKIKQLLQENEV QFLDKGDENTVIDLGSKRQY ADCSEIFNDGYKLSGFYKIK PLQSPAEFVYCDMSDGGGW TVIQRSDGSENFNRGWKDY ENGFGNFVQKHGEYWLGNKN LHFLTQTQEDYTLKIDLADFE KNSRYAQYKNFKVGDEKNFY ELNIGEYSGTAGDSLGNFH PEVQWWASHQRMKFSTWDR