

Anti-UGP2 antibody (20-250 aa) (STJ11107404)

STJ11107404

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:200-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	7360
Gene Symbol	UGP2
Uniprot ID	UGPA_HUMAN
Immunogen	
Immunogen Region	20-250 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 20-250 of human UGP2 (NP_006750.3).
Immunogen Sequence	FQEVIRQEELSVKKELEKI LTTASSHEFEHTKKDLGFR KLFHRFLQEKGPSVDWGKIQ RPPEDSIQPYEKIKARGLPD NISSVLNKLNVVVLNGGLGT SMGCKGPKSLIGVRNENTFL DLTVQQIEHLNKTYNQDVPL VLMNSFNTDEDTKKILQKYN HCRVKIYTFNQSRYPRIKE SLLPVAKDVSYSGENTEAWY PPGHGDYASFYNSGLLDTF IGEKKEYIFVS