

Anti-AKAP10 antibody (413-662 aa) (STJ11107407)

STJ11107407

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	11216
Gene Symbol	AKAP10
Uniprot ID	AKA10_HUMAN
Immunogen	
Immunogen Region	413-662 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 413-662 of human AKAP10 (NP_009133.2).
Immunogen Sequence	DNFQSQLAAKKGQYDGQEAQ NDAMILYDKYFSLQATHPLG FDDVVRLEIESNICREGGPL PNCFTTPLRQAWTTMEKVFL PGFLSSNLYYKYLNDLIHSV RGDEFLGGNVSLTAPGSVGP PDESHPGSSDSSASQSSVKK ASIKILKNFDEAIIVDAASL DPESLYQRTYAGKMTFGRVS DLGQFIRESEPEPDVRKSKG SMFSQAMKKWVQGNTDEAQE ELAWKIAKMIVSDIMQQAQ