

Anti-HSPBP1 antibody (1-359 aa) (STJ11107408)

STJ11107408

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:200-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	23640
Gene Symbol	HSPBP1
Uniprot ID	HPBP1_HUMAN
Immunogen	
Immunogen Region	1-359 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-359 of human HSPBP1 (NP_036399.3).
Immunogen Sequence	MSDEGSRGSRLPLALPPASQ GCSSGGGGGGSSAGGSGNSR PPRNLQGLLQMAITAGSEEP DPPPEPMSEERRQWLQEAMS AAFRGQREEVEQMKSLRVL SQMPMPPTAGEAEQAADQQR EGALELLADLCENMDNAADF CQLSGMHLLVGRYLEAGAAG LRWRAAQLIGTCSQNVAAIQ EQVLGLGALRKLLRLDRDA CDTVVRKALFAISCLVREQE AGLLQFLRLDGFSLMRAM