

Anti-CNNM3 antibody (410-640 aa) (STJ11107409)

STJ11107409

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	26505
Gene Symbol	CNNM3
Uniprot ID	CNNM3_HUMAN
Immunogen	
Immunogen Region	410-640 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 410-640 of human CNNM3 (NP_951060.1).
Immunogen Sequence	DTVVKRKPPASLMAPLKREE FSLFKVSDDEYKVTISPQLL LATQRFLSREVDVFSPRLIS EKVLLHLLKHPSVNQEVRFD ESNRLATHHYLYQRSQPVVDY FILILQGRVEIGKEGLKF ENGAFTYYGVSALTPSSVH QSPVSSLQPIRHDLQPDPGD GTHSSAYCPDYTVRALSDLQ LIKVTRLQYLNALLATRAQN LPQSPENTDLQVIPGSQTRL LGEKTTTAAGS