

Anti-BBS9 antibody (55-300 aa) (STJ11107410)

STJ11107410

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:200-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	27241
Gene Symbol	BBS9
Uniprot ID	PTHB1_HUMAN
Immunogen	
Immunogen Region	55-300 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 55-300 of human BBS9 (NP_055266.2).
Immunogen Sequence	HPAKTGDGAQAEDLLLEVDL RDPVLQVEVGKFSVSGTEMLH LAVLHSRKL CVYSVSGTLGN VEHGNQCQMKLMYEHNLRQT ACNMTYGSFGGVKGRDLICI QSMGMLMVFEQESYAFGRF LPGFLLPGPLAYSSRTDSFL TVSSCQQVESYKYQVLAFAT DADKRQETEQKLGSGKRLV VDWTNLNIGEQALDICVSFN QSASSVFLGERNFFCLKDN GQIRFMKKLDWSPSCFLPY