

Anti-FAM175A antibody (1-300 aa) (STJ11107413)

STJ11107413

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	84142
Gene Symbol	ABRAXAS1
Uniprot ID	ABRX1_HUMAN
Immunogen	
Immunogen Region	1-300 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human FAM175A (NP_620775.2).
Immunogen Sequence	MEGESTSAVLSGFVLGALAF QHLNTDSDTEGFLGGEVKGGE AKNSITDSQMDDVEVYTID IQKYIPCYQLFSFYNSSGEV NEQALKKILSNVKKNVVGWY KFRRHSDQIMTFRERLLHKN LQEHFSNQDLVFLLLTPSII TESCSTHRLHSLYKPKQKGL FHRVPLVANLGMSEQLGYK TVSGSCMSTGFSRAVQTHSS KFFEEDGSLKEVHKINEMYA SLQEELKSICKKVEDSEQA