

Anti-SH2D1B antibody (1-132 aa) (STJ11107414)

STJ11107414

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:200-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	117157
Gene Symbol	SH2D1B
Uniprot ID	SH21B_HUMAN
Immunogen	
Immunogen Region	1-132 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-132 of human SH2D1B (NP_444512.2).
Immunogen Sequence	MDLPYYHGRLTKQDCETLLL KEGVDGNFLLRDSEIPGVL CLCVSFKNIVYTYRIFREKH GYYRIQTAEGSPKQVFPCLK ELISKFEKPNQGMVVHLLKP IKRTSPSLRWRLKLELETF VNSNSDYVDVLP