

Anti-SLC12A1 antibody (1-177 aa) (STJ11107416)

STJ11107416

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	6557
Gene Symbol	SLC12A1
Uniprot ID	S12A1_HUMAN
Immunogen	
Immunogen Region	1-177 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-177 of human SLC12A1 (NP_000329.2).
Immunogen Sequence	MSLNSSNVFLDSVPSNTNR FQVSVINENHESSAAADDNT DPPHYEETSFGDEAQKRLRI SFRPGNQECYDNFLQSGETA KTDASFHAYDSHTNTYYLQT FGHNTMDAVPKIEYYRNTGS ISGPKVNRPSLLEIHEQLAK NVAVTPSSADRVANGDGIPG DEQAENKEDDQAGVVKF