

Anti-KIAA0020 antibody (1-220 aa) (STJ11107423)

STJ11107423

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	9933
Gene Symbol	PUM3
Uniprot ID	PUM3_HUMAN
Immunogen	
Immunogen Region	1-220 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human KIAA0020 (NP_055693.4).
Immunogen Sequence	MEVKGKKQFTGKSTKTAQEK NRFHKNSDSGSSKTFPTRKV AKEGGPKVTSRNFESITKL GKKGVKQFKNKQQGDKSPKN KFQPANFNNKRFQPDGRS DESAAKPKWDDFKKKKKEL KQSRQLSDKTNVDIVRAKQ MWEILRRKDCDKEKRVKLMS DLQKLIQGIKTIAFAHDST RVIQCYIQYGNEEQRKQAFE ELRDDLVLSKAKYSRNIVK