

Anti-KDM3B/KDM3A antibody (1472-1761 aa) (STJ11107426)

STJ11107426

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	55818 51780
Gene Symbol	KDM3A KDM3B
Uniprot ID	KDM3A_HUMAN KDM3B_HUMAN
Immunogen	
Immunogen Region	1472-1761 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1472-1761 of human KDM3B (NP_057688.2).
Immunogen Sequence	RLRSEDGQPMVLKLDWPPG EDFRDMMPTRFEDLMENLPL PEYTKRDGRLNLASRLPSYF VRPDLGPKMYNAYGLITAED RRVGTTNLHLDVSDAVNVMV YVGIPIGEGAHDEEVLTID EGDDEVTKQRIHDGKEKPG ALWHIYAAKDAEKIRELLRK VGEEQGQENPPDHPIDQS WYLDQTLRKRLYEYGVQGW AIVQFLGDAVFIPAGAPHQV HNLYSCIKVAEDFVSPEHV