

Anti-CIDECD antibody (1-238 aa) (STJ11107429)

STJ11107429

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	63924
Gene Symbol	CIDECD
Uniprot ID	CIDECD_HUMAN
Immunogen	
Immunogen Region	1-238 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-238 of human CIDECD (NP_001186481.1).
Immunogen Sequence	MEYAMKSLSLYPKSLSRHV SVRTSVVTQQLSEPSPKAP RARPCRVSTADRSVRKGIMA YSLEDLLLVKVRDTLMLADKP FFLVLEEDGTTVETEEYFQA LAGDTVMVLQKGQKWQPPS EQGTRHPLSLSHKPAKKIDV ARVTFDLYKLNPDQDFIGCLN VKATFYDTYSLSYDLHCCGA KRIMKEAFRWALFSMQATGH VLLGTSCYLQQLDATEEGQ PPKGKASSLIPTCLKILQ