

Anti-TSPAN32 antibody (110-205 aa) (STJ11107433)

STJ11107433

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	10077
Gene Symbol	TSPAN32
Uniprot ID	TSN32_HUMAN
Immunogen	
Immunogen Region	110-205 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 110-205 of human TSPAN32 (NP_620591.3).
Immunogen Sequence	FWRLHSPTQVEDAMLDTYDL VYEQAMKGTSHVRRQELAAI QDVFLCCGKKSPFSRLGSTE ADLCQGEEAAREDCLQGIRS FLRTHQQVASSLTSG

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081