

Anti-UBE3C antibody (970-1083 aa) (STJ11107434)

STJ11107434

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IF/ICC:1:10-1:100 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	9690
Gene Symbol	UBE3C
Uniprot ID	UBE3C_HUMAN
Immunogen	
Immunogen Region	970-1083 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 970-1083 of human UBE3C (NP_055486.2).
Immunogen Sequence	PISLEDLKSFTNYSGGYSAD HPVVKVFWRVVEGFTDEEKR KLLKFVTSCSRPPLLGFKEL YPAFCIHNGGSDLERLPTAS TCMNLKLPFYDETLLRSK LLYAIECAAGFELS