

Anti-FRS2 antibody (200-470 aa) (STJ11107437)

STJ11107437

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	10818
Gene Symbol	FRS2
Uniprot ID	FRS2_HUMAN
Immunogen	
Immunogen Region	200-470 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 200-470 of human FRS2 (NP_006645.3).
Immunogen Sequence	TGVQEEERKNRTSVHVPLEAR VSNAESSTPKKEEPSSIEDRD PQILLEPEGVKFVLGPTPVQ KQLMEKEKLEQLGRDQVSGS GANNTEWDTGYDSDEERRDAP SVNKLVEYENINGLSIPSASG VRRGRLTSTSTSDTQNNNS AQRRTALLNYENLPSLPPVW EARKLSRDEDDNLGPKTPSL NGYHNNLDPMHNYVNTENVT VPASAHKIEYSRRRDCTPTV FNFDIRRPSLEHRQLNYIQ