

**Anti-GLRX2 antibody (1-124 aa) (STJ11107440)**

STJ11107440

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	
<b>Applications</b>	WB/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	Lot specific
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
<b>Formulation</b>	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">51022</a>
<b>Gene Symbol</b>	<a href="#">GLRX2</a>
<b>Uniprot ID</b>	<a href="#">GLRX2_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	1-124 aa
<b>Specificity</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 1-124 of human GLRX2 (NP_001230328.1).
<b>Immunogen Sequence</b>	MESNTSSLENLATAPVNQI QETISDNCVVIFSKTSCSYC TMAKKLFHDMNVNYKVVELD LLEYGNQFQDALYKMTGERT VPRIFVNGTFIGGATDTHRL HKEGKLLPLVHQCYLKKSQR KEFQ