

Anti-GNLY antibody (23-145 aa) (STJ11107441)

STJ11107441

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	10578
Gene Symbol	GNLY
Uniprot ID	GNLY_HUMAN
Immunogen	
Immunogen Region	23-145 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 23-145 of human GNLY (NP_006424.2).
Immunogen Sequence	RLSPEYYDLARAHLRDEEKS CPCLAQEGPQGDLLTKTQEL GRDYRTCLTIVQKLKKMVDK PTQRSVSNAATRVCGRTGRSR WRDVCRNFMRRYQSRVTQGL VAGETAQQICEDLRLCIPST GPL