

Anti-PPP4R2 antibody (188-417 aa) (STJ11107447)

STJ11107447

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	151987
Gene Symbol	PPP4R2
Uniprot ID	PP4R2_HUMAN
Immunogen	
Immunogen Region	188-417 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 188-417 of human PPP4R2 (NP_777567.1).
Immunogen Sequence	MTTNGLPESDTSKEANLQQN EEKNHSDSSTSESEVSSVSP LKNKHPDEDAVEAGHEVKR LRFDKGEVRETASQTTSSSE ISSVMVGETEASSSSQDKDK DSRCTRQHCTEDEEEDEEEE EEESFMTSREMIPERKNQEK ESDDALTVNEETSEENNQME ESDVSQAEKDLLHSEGENE GPVSSSSSDCRETELVGSN SSKTGEILSESSMENDDEAT EVTDEPMEQD