

## Anti-CSN1S1 antibody (16-135 aa) (STJ11107451)

STJ11107451

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	
<b>Applications</b>	WB/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	Lot specific
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
<b>Formulation</b>	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	1446
<b>Gene Symbol</b>	CSN1S1
<b>Uniprot ID</b>	CASA1_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	16-135 aa
<b>Specificity</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 16-135 of human CSN1S1 (NP_001881.1).
<b>Immunogen Sequence</b>	RPKLPLRYPERLQNPSESSE PIPLESREEYMNGMNRQRNI LREKQTDEIKDTRNESTQNC VVAEPEKMESSISSSEEMS LSKCAEQFCRLNEYNQLQLQ AAHAQEQIRRMNENSHVQVP