

Anti-Mammaglobin A antibody (1-93 aa) (STJ11107452)

STJ11107452

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	4250
Gene Symbol	SCGB2A2
Uniprot ID	SG2A2_HUMAN
Immunogen	
Immunogen Region	1-93 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-93 of human Mammaglobin A (NP_002402.1).
Immunogen Sequence	MKLLMVLMLAALSQHCYAGS GCPLLENVISKTINPQVSKT EYKELLQEFIDNATTNAID ELKECFLNQTDETLSNVEVF MQLIYDSSLCDLF