

Anti-PRND antibody (20-160 aa) (STJ11107454)

STJ11107454

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	23627
Gene Symbol	PRND
Uniprot ID	PRND_HUMAN
Immunogen	
Immunogen Region	20-160 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 20-160 of human PRND (NP_036541.2).
Immunogen Sequence	HL SAVQTRGIKHRIKWNRKA LPSTAQITEAQVAENRPGAF IKQGRKLDIDFGAEGNRYYE ANYWQFPDGIHYNGCSEANV TKEAFVTGCINATQAANQGE FQKPDNKLHQQVLWRLVQEL CSLKHCFWLERGAGLRVTM H