

Anti-CDR2 antibody (175-454 aa) (STJ11107457)

STJ11107457

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	1039
Gene Symbol	CDR2
Uniprot ID	CDR2_HUMAN
Immunogen	
Immunogen Region	175-454 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 175-454 of human CDR2 (NP_001793.1).
Immunogen Sequence	HVFAEKITSLQGQSPDEEE NEHLKKTVTMLQAQLSLERQ KRVTMEEYGLVLKENSELE QQLGATGAYRARALELEAEV AEMRQMLQSEHPFVNGVEKL VPDSLYVPFKEPSQSLLEEM FLTVPESHKPLKRSSSETI LSSLAGSDIVKGHEETCIRR AKAVKQRGISLLHEVDQYS ALKVKEYELLKKCQEEQDSL SHKAVQTSRAAKDLTGVNA QSEPVASGWELASVNPPEV