

Anti-CSRNP2 antibody (324-543 aa) (STJ11107458)

STJ11107458

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	81566
Gene Symbol	CSRNP2
Uniprot ID	CSRNP2_HUMAN
Immunogen	
Immunogen Region	324-543 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 324-543 of human CSRNP2 (NP_110436.1).
Immunogen Sequence	VMHLQSAEELERLKAEDSS GSSASLDSSIESLGVCIIEE PLAVPEELCPGLTAPILIA QLPPGSSVLCFTENSDHPTA STVNSPSYLSNGPLVYYQVE QRPVLGVKGEPGTEEGSASF PKEKDLNVFSLPVTSLVACS STDPAALCKSEVGKTPITLEA LLPEDCNPEEPENEDFHPSW SPSSLPFRTDNEEGCGMVKT SQQNEDRPPEDSSLELPLAV