

Anti-WDR47 antibody (350-600 aa) (STJ11107464)

STJ11107464

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	22911
Gene Symbol	WDR47
Uniprot ID	WDR47_HUMAN
Immunogen	
Immunogen Region	350-600 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 350-600 of human WDR47 (NP_001136023.1).
Immunogen Sequence	YPGVQNLSRSLMLETECHS IYEESEPERDTPVDAQRPIS EILGQSSVSEKEPANGAQNPGPAKQEKNELRDSTEQFQY YRQRLRYQQHLEQKEQQRQI YQQLMLLEGGVNQEDGPDQQQ NLTEQFLNRSIQKLGLNIG MDGLGNEVSALNQQCNGSKG NGSNGSSVTSFTTPPDSSQ RLTHDASNIHTSTPRNPGST NHIPFLEESPCGSQISSEHS VIKPPLGDSPGSLRSKGE