

Anti-PLEKHF2 antibody (1-249 aa) (STJ11107467)

STJ11107467

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	79666
Gene Symbol	PLEKHF2
Uniprot ID	PKHF2_HUMAN
Immunogen	
Immunogen Region	1-249 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-249 of human PLEKHF2 (NP_078889.1).
Immunogen Sequence	MVDRLANSEANTRRISIVEN CFGAAGQPLTIPGRVLIGEG VLTCLCRKKPKARQFFLFND ILVYGNIVIQKKKYNKQHII PLENVTIDSIKDEGDLRNGW LIKTPTKSFAYAAATEKS EWMNHINKCVTDLLSKSGKT PSNEHAAVWVPDSEATVCMR CQKAKFTPVNRRHHCRKCGF VCGPCSEKRFLPSQSSKP VRICDFCYDLLSAGDMATCQ PARSDSYSQSLKSPLNDMS