

Anti-ERMP1 antibody (673-904 aa) (STJ11107468)

STJ11107468

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	79956
Gene Symbol	ERMP1
Uniprot ID	ERMP1_HUMAN
Immunogen	
Immunogen Region	673-904 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 673-904 of human ERMP1 (NP_079172.2).
Immunogen Sequence	PYSSNPANPKPKRVFLQHMT RTFHDLEGNAVKRDSGIWIN GFDYTGISHITPHIPEINDS IRAHCEENAPLCGFPWYLPV HFLIRKNWYLPAPVSPRNP PHFRLISKEQTPWDSIKLTF EATGPSHMSFYVRAHKGSTL SQWSLGNTPVTSTKGGDYFV FYSHGLQASAWQFWIEVQVS EEHPEGMTVAIAAHYLSGE DKRSPQLDALKEKFPDWTFP SAWVCTYDLFVF