

Anti-AMY2B antibody (16-320 aa) (STJ11107471)

STJ11107471

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	280
Gene Symbol	AMY2B
Uniprot ID	AMY2B_HUMAN
Immunogen	
Immunogen Region	16-320 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 16-320 of human AMY2B (NP_066188.1).
Immunogen Sequence	QYSPNTQQGRTSIVHLFEWR WVDIAECERYLAPKGFQGV QVSPPNENVIAHNPFPPWWE RYQPVSYKLCTRSGNEDEFR NMVTRCINNNGVRIYDAVIN HMSGNAVSAGTSSTCGSYFN PGSRDFPAVPYSGWDFNDGK CKTGSGDIENYNDATQVRDC RLVGLLDLALEKDYVRKIA EYMNHLIDIGVAGFRLDASK HMWPGDIKAILDKLHNLNSN WFPAGSKPFIYQEVIDLGG