

Anti-MOCS1 antibody (90-360 aa) (STJ11107473)

STJ11107473

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	4337
Gene Symbol	MOCS1
Uniprot ID	MOCS1_HUMAN
Immunogen	
Immunogen Region	90-360 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 90-360 of human MOCS1 (NP_005934.2).
Immunogen Sequence	EEGVPLTPKANLLTTEEILT LARLFVKEGIDKIRLTGGEP LIRPDVVDIVAQLQRLEGLR TIGVTTNGINLARLLPQLQK AGLSAINISLDTLVPKFEF IVRRKGFHKVMEGIHKAIEL GYNPVKVNCCVMRGLNEDEL LDFAALTEGLPLDVRFIEYM PFDGNKWNFKKMVSYPEMLD TVRQQWPELEKVPPEESSTA KAFKIPGFQGGISFITS MSE HFCGTCNRLRITADGNLKV