

Anti-PRDM10 antibody (600-870 aa) (STJ11107475)

STJ11107475

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	56980
Gene Symbol	PRDM10
Uniprot ID	PRD10_HUMAN
Immunogen	
Immunogen Region	600-870 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 600-870 of human PRDM10 (NP_064613.2).
Immunogen Sequence	LDLLKDHVAIHINDGYFTCP TCKKRFPDFIQVKKHVRSFH SEKIYQCTECDKAFCRPDKL RLHMLRHSDRKDFLCSTCGK QFKRKDKLREHMQRMHNP ER EAKKADRISRSTFKPRITS TDYDSFTFKCRLCMMGFRRR GMLVNHLSKRHPDMKIEEVP ELTLPIIKPNRDYFCQYCDK VYKSASKRKAHILKNHPGAE LPPSIRKLRPAGPGEPDPM L STHTQLTGTIATPPVCCPH