

Anti-FASTKD3 antibody (1-220 aa) (STJ11107476)

STJ11107476

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	79072
Gene Symbol	FASTKD3
Uniprot ID	FAKD3_HUMAN
Immunogen	
Immunogen Region	1-220 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human FASTKD3 (NP_076996.2).
Immunogen Sequence	MALITLRKNLYRLSDFQMHR ALAALKNKPLNHVHKVVKER LCPWLCSRQPEPFGVKFHHH HCKKFHSHKNGNDLHPLGGPV FSQVSDCDRLQNVKNEESQ MFYRRLSNLTSSEEVLSFIS TMTLPDPTMAAGALQRICEV EKKDGDQGLPKEILENSIFQ ALCFQFEKEPSQLSNTSLVT ALQALILLHVDPPQSLLLLNL VAECQNRLRKGGMEVRNLGI