

Anti-PHEX antibody (520-749 aa) (STJ11107479)

STJ11107479

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	5251
Gene Symbol	PHEX
Uniprot ID	PHEX_HUMAN
Immunogen	
Immunogen Region	520-749 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 520-749 of human PHEX (NP_000435.3).
Immunogen Sequence	WLRKAVPKTEWFTNPPTVNA FYSASTNQIRFPAGELQKPF FWGTEYPRSLSYGAIGVIVG HEFTHGFDNNGRKYDKNGNL DPWWSTESEEFKEKTKCMI NQYSNYYWKKAGLNVKGKRT LGENIADNGGLREAFRAYRK WINDRRQGLEELPLPGITFT NNQLFFLSYAHVRCNSYRPE AAREQVQIGAHSPQFRVNG AISNFEEFQKAFNCPNSTM NRGMDSCRLW