

Anti-FOXP3 antibody (1-100 aa) (STJ11107482)

STJ11107482

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION


Gene ID	50943
Gene Symbol	FOXP3
Uniprot ID	FOXP3_HUMAN
Immunogen	
Immunogen Region	1-100 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FOXP3 (NP_054728.2).
Immunogen Sequence	MPNPRPGKPSAPSLALGPSP GASPSWRAAPKASDLLGARG PGGTFQGRDLRGGAHASSSS LNPMPPSQLQLPTLPLVMVA PSGARLGLPLHLQALLQDRP



Immunohistochemistry analysis of paraffin-embedded Human tonsil using FOXP3 Rabbit pAb (STJ11107482) at dilution of 1:100 (40x lens). Microwave antigen retrieval performed with 0.01M PBS Buffer (pH 7.2) prior to immunohistochemistry staining.



Immunohistochemistry analysis of paraffin-embedded Human appendix using FOXP3 Rabbit pAb (STJ11107482) at dilution of 1:100 (40x lens). Microwave antigen retrieval performed with 0.01M PBS Buffer (pH 7.2) prior to immunohistochemistry staining.



Immunohistochemistry analysis of paraffin-embedded Human stomach cancer using FOXP3 Rabbit pAb (STJ11107482) at dilution of 1:100 (40x lens). Microwave antigen retrieval performed with 0.01M PBS Buffer (pH 7.2) prior to immunohistochemistry staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081