

Anti-EDN2 antibody (69-178 aa) (STJ11107483)

STJ11107483

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	1907
Gene Symbol	EDN2
Uniprot ID	EDN2_HUMAN
Immunogen	
Immunogen Region	69-178 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 69-178 of human EDN2 (NP_001947.1).
Immunogen Sequence	WVNTPEQTAPYGLGNPPRRR RRS�PRRCQCSARDPACAT FCLRRPWTEAGAVPSRKSPA DVFQTGKTGATTGELLQRLR DISTVKSLFAKRQQEAMREP RSTHSRWKR