

Anti-AIMP2 antibody (141-320 aa) (STJ11107491)

STJ11107491

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	7965
Gene Symbol	AIMP2
Uniprot ID	AIMP2_HUMAN
Immunogen	
Immunogen Region	141-320 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 141-320 of human AIMP2 (NP_006294.2).
Immunogen Sequence	LLCEHFRVLSTVHTHSSVKS VPENLLKCFGEQNKQPRQD YQLGFTLIWKNVPKTQMKFS IQTMCPIEGEGNIARFLFSL FGQKHNAVNATLIDSWVDIA IFQLKEGSSKEAAVFRSMN SALGKSPWLAGNELTVADW LWSVLQIGGCSVTVPANVQ RWMRSCENLAPFNTALKLLK