

Anti-COPS4 antibody (1-310 aa) (STJ11107498)

STJ11107498

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	51138
Gene Symbol	COPS4
Uniprot ID	CSN4_HUMAN
Immunogen	
Immunogen Region	1-310 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-310 of human COPS4 (NP_057213.2).
Immunogen Sequence	MAAAVRQDLAQLMNSSGSHK DLAGKYRQILEKAIQLSGAE QLEALKAFVEAMVNVNSLV ISRQLLTDFCTHLPNLPDST AKEIYHFTLEKIQPRVISFE EQVASIRQHLASIYEKEEDW RNAAQVLVGIPLETGQKQYN VDYKLETYLKARLYLEDD PVQAEAYINRASLLQNESTN EQLQIHYKVCYARVLDYRRK FIEAAQRYNELSYKTIVHES ERLEALKHALHCTILASAG