

Anti-Renalase antibody (31-342 aa) (STJ11107501)

STJ11107501

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	55328
Gene Symbol	RNLS
Uniprot ID	RNLS_HUMAN
Immunogen	
Immunogen Region	31-342 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 31-342 of human Renalase (NP_001026879.2).
Immunogen Sequence	AVWDKAEDSGGRMTTACSPH NPQCTADLGAQYITCTPHYA KKHQRFYDELLAYGVLRLPLS SPIEGMVMKEGDCNFBAPQG ISSIIKHYLKESGAEVYFRH RVTQINLRDDKWEVSKQTGS PEQFDLVLTMPVPEILQLQ GDITTLISECQRQQLAVSY SSRYALGLFYEAGTKIDVPW AGQYITSNPCIRFVSIDNKK RNIESSEIGPSLVIHTTVPF GVTYLEHSIEDVQELVFQQ