

Anti-IL22RA1 antibody (16-230 aa) (STJ11107505)

STJ11107505

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	58985
Gene Symbol	IL22RA1
Uniprot ID	I22R1_HUMAN
Immunogen	
Immunogen Region	16-230 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 16-230 of human IL22RA1 (NP_067081.2).
Immunogen Sequence	HAPEDPSDLLQHVKFQSSNF ENILTWDSGPEGTPDVTYSI EYKTYGERDWWAKKGCCRIT RKSCNLTVETGNLTELYAR VTAVSAGGRSATKMTDRFSS LQHTTLKPPDVTCSISKVRSI QMIVHPTPTPIRAGDGHRLT LEDIFHDLFYHLELQVNRTY QMHLGGKQREYEFFGLTPDT EFLGTIMICVPTWAKESAPY MCRVKTLPDRTWTYS