

Anti-EDA2R antibody (1-138 aa) (STJ11107507)

STJ11107507

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	60401
Gene Symbol	EDA2R
Uniprot ID	TNR27_HUMAN
Immunogen	
Immunogen Region	1-138 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-138 of human EDA2R (NP_068555.1).
Immunogen Sequence	MDCQENEYWDQWGRVCVTCQR CGPGQELSKDCGYGEGGDY CTACPPRRYKSSWGHHCQS CITCAVINRVQKVNCTATSN AVCGDCLPRFYRKTRIGGLQ DQECIPCTKQTPTEVQCAF QLSLVEADPTVPPQEAT