

Anti-SPRED2 antibody (120-300 aa) (STJ11107514)

STJ11107514

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	200734
Gene Symbol	SPRED2
Uniprot ID	SPRE2_HUMAN
Immunogen	
Immunogen Region	120-300 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 120-300 of human SPRED2 (NP_861449.2).
Immunogen Sequence	EDLIEGSTTSSSTIHNEAL GDDDVFTTATDSSSNSSQKR EQPTRTISSPTSCEHRRITYT LGHLHDSYPTDHYHLDQPMR RPYRQVSFPDDDEIVRINP REKIWMGTGYEDYRHAPVRGK YDPSEDADSSYVRFKGEV PKHDYNYPYVDSSDFGLGED PKGRGGSVIKTQPSRGKSRR R