

Anti-GABRA1 antibody (28-251 aa) (STJ11107518)

STJ11107518

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	2554
Gene Symbol	GABRA1
Uniprot ID	GBRA1_HUMAN
Immunogen	
Immunogen Region	28-251 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 28-251 of human GABRA1 (NP_000797.2).
Immunogen Sequence	QPSLQDELKDNTTVFTRILD RLLDGYDNLRPLGERVTE VKTDIFVTSFGPVSDHMEY TIDVFFRQSWKDERLKFKGP MTVLRLLNNLMASKIWTPDTF FHNGKKSVAHNMTMPNKLRL ITEDGTLTYTMRLTVRAECP MHLEDFPMDAHACPLKFGSY AYTRAENVYEWTPARSVV VAEDGSRLNQYDLLGQTVDS GIVQSSTGEYVWMTTHFLK RKIG