

Anti-MEX3C antibody (484-659 aa) (STJ11107520)
STJ11107520

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	51320
Gene Symbol	MEX3C
Uniprot ID	MEX3C_HUMAN
Immunogen	
Immunogen Region	484-659 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 484-659 of human MEX3C (NP_057710.3). LPSVGSSEDLAVDSPAFDSLPTSAQTIWTPFEPVNPVLSGFGSDPSGNMKTQRRGSQPSLTPRLSPTFPESIEHPLARRVRSD
Immunogen Sequence	PPSTGNHVLPIYIPAFSNGTNSYSSSNGGSTSSPPESRKHDCVICFENEVIAALVPCGHNLCFMECANKICEKRTPS CPVCQTAVTQAIQIHS



Western blot analysis of lysates from SH-SY5Y cells, using MEX3C Rabbit pAb (STJ11107520) at 1:1000 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution.
Lysates/proteins: 25 Mu g per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit
Exposure time: 15s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081