

Anti-BRF1 antibody (1-208 aa) (STJ11107523)

STJ11107523

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	2972
Gene Symbol	BRF1
Uniprot ID	TF3B_HUMAN
Immunogen	
Immunogen Region	1-208 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-208 of human BRF1 (NP_001229719.1).
Immunogen Sequence	MTGRVCRGCGGTDIELDAAR GDAVCTACGSVLEDNIIVSE VQFVSSGGGSSAVGQFVSL DGAGKTPTLGGGFHVNLGKE SRAQTLQNGRRHHHLGNQL QLNQHCLDTAFNFFKMAVSR HLTRGRKMAHVIAACLYLVC RTEGTPHMLLDLSDLLQVDS LRPASFTWGC DLGVVTRVV TGVYPRCLHASQWPVCAACP VRKFWSVG