

Anti-JCHAIN antibody (23-159 aa) (STJ11107532)

STJ11107532

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	3512
Gene Symbol	JCHAIN
Uniprot ID	IGJ_HUMAN
Immunogen	
Immunogen Region	23-159 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 23-159 of human JCHAIN (NP_653247.1).
Immunogen Sequence	QEDERIVLDNKCKCARITS RIIRSSDPNEDIVERNIRI IVPLNNRENISDPTSPLRTR FVYHLSDLCKKCDPTEVELD NQIVTATQSNICDEDSATET CYTYDRNKCYTAVVPLVYGG ETKMVETALTPDACYPD