

Anti-SLCO2A1 antibody (420-520 aa) (STJ11107533)

STJ11107533

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	6578
Gene Symbol	SLCO2A1
Uniprot ID	SO2A1_HUMAN
Immunogen	
Immunogen	420-520 aa
Region	
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 420-520 of human SLCO2A1 (NP_005621.2).
Immunogen	CSTPTVAEVYPPSTSSSIHP QSPACRRDCSPDSIFHPVC GDNGIEYLSPOCHAGCSNINM SSATSKQLIYLNCSCVTGGS
Sequence	ASAKTGSCPVPCAHLFLPAI F