

Anti-TBCD antibody (1-240 aa) (STJ11107535)

STJ11107535

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:200-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	6904
Gene Symbol	TBCD
Uniprot ID	TBCD_HUMAN
Immunogen Region	1-240 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human TBCD (NP_005984.3).
Immunogen Sequence	MALSDEPAAGGPEEEAEDET LAFGAALFAFGESAETRAL GRLREVHGGGAEREVALERF RVIMDKYQEQPHLLDPHLEW MMNLLLIVQDQTSPASLVH LAFKFLYIITKVRGYKTF LR LFPHEVADVEPVLDLVTIQN PKDHEAWETRYMLLLWLSVT CLIPDFSRDLGNLLTQPGQ ARMSIMDRILQIAESYLIVS DKARDAAAVLVSRFITRPDV KQSKMAEFLDWSLCNLARS